

Housewares & Furnishings

PikStik Pro

- Lightweight design is essential for those with limited strength.
- No sharp edges or ways to pinch fingers means the product is usable by all, even kids.
- Whether at home, work, in the yard, or in water - this design will last and not rust.
- Free replacement within one year of purchase, if damage is a result of normal wear and tear.
- Made Out of GE's "Lexan" Material



Yellow

20" - SKU: PS-201

26" - SKU: PS-261

32" - SKU: PS-321

48" - SKU: PS-488

Orange

20" - SKU: PS-204

26" - SKU: PS-264

32" - SKU: PS-324

48" - SKU: PS-484



32" SKU: PS-L320

PikStik Lite

- Jaws can be locked in any position from the handle without having to reach a cumbersome locking mechanism with your thumb. You no longer have to grip and then lock, just keep gripping for an automatic lock!
- Jaws rotate 360 degrees; reach and pick up objects from any angle.
- Dual hinged lever creates greater jaw gripping power without needing extra hand strength.
- Four finger rests make it easier and more comfortable to grip.

PikStik Industrial

- Wide Jaw Opening: Pick up more varied and larger objects.
- Tips will never fall off, they are molded right onto the jaws.
- Unique dual vault channel design adds unprecedented strength to shaft.
- Locking mechanism at the base of the handle gives you the leverage of your entire hand to lock the jaws, not just your thumb! This makes it easily usable by all (right or left handed), even those without a lot of strength in their hands.



Yellow

32" - SKU: PS-I-321

48" - SKU: PS-I-481

email: sales@egeneralmedical.com

phone: 919.844.9402

fax: 919.844.9403

Raleigh, NC 27616 USA

