

# HappyLight® | Deluxe

Sunshine Supplement Light System



- **Brighter Natural Spectrum® Light:** Dual patented 36 W high-output, 5,000-hour bulbs with Tricolite Phosphor Technology™ deliver therapeutic 10,000 LUX of light without glare.
- **Versatile Design:** Weighs just 7 lbs. and comes with built in handle and table-top stand plus space-saving wall-mount option.
- **Best-in-class Power System:** Electronic Ballast and Flicker Elimination Technology™ assure no-buzz operation and instant-on, flicker-free illumination without generating EMFs.
- **Guaranteed Performance:** Satisfaction guaranteed by Verilux® — for 50 years the leaders in Natural Spectrum® indoor sunlight — and backed by 1-year limited warranty.




# HappyLight® | Deluxe

## Sunshine Supplement Light System

10,000 LUX brightens your life! A must for everyone who spends a lot of time indoors, this portable unit simulates sunshine to activate your own natural mood enhancers without exposing you to harmful UV rays or EMFs.

A patented Verilux® bulb delivers the 10,000 LUX of healthful light the National Institute of Mental Health recommends for treating Seasonal Affective Disorder (S.A.D.). HappyLight® provides this therapy-strength Natural Spectrum® light using exclusive Verilux® Flicker-Elimination Technology™.

Unit can be wall-mounted or set on a tabletop so it's easy to get the nourishing light you need for brighter moods, greater well-being, and clearer thinking.



## Specifications

### Models - Colors

#### HPLD

Lens	Frosted
Body	Ivory
Cord/Plug	White
UPC	7 68533 01036 0

#### Ballast

Ballast Type	High Output Electronic
--------------	------------------------

#### Packaging

##### Unit Carton Box (UCB)

Quantity	1
Dimensions	22"D x 16.125"H x 5.125"W
Gross Weight (GW)	8.4 lbs

##### Outer Unit Shipper (OUS)

Quantity	1
Dimensions	23" x 17" x 5.625"
Gross Weight (GW)	9.7 lbs

### Product

#### Overall Dimensions

Weight	6.8 lbs
Dimensions:	19.5"H x 13.375"W x 3.25"D

#### Other

Body Material	ABS
Reflector	Specular Metal

#### Switch

Switch Type	On/Off
Switch Location	Top of Body

#### Cord/Plug

Plug Type	USA three prong (grounded)
Cord Length	96" (8 feet)
Proposition 65	Compliant

#### Bulb

CFL36VLX	36 Watt Natural Daylight Fluorescent
----------	--------------------------------------

#### Safety

UL Listed	
-----------	--